

Model 14014

STAND PERIMETER SCHWEIGGER TYPE USER INSTRUCTIONS



3700 Sagamore Parkway North
P.O. Box 5729 • Lafayette, IN 47903 USA
Tel: 765.423.1505 • 800.428.7545
Fax: 765.423.4111
E-mail: info@lafayetteinstrument.com
www.lafayetteinstrument.com

S.I. Instruments
256 South Rd. Hilton
South Australia 5033
Ph (08) 8352 5511

info@si-instruments.com
www.si-instruments.com

Table of Contents:

System Description	3
Procedure	3
Table 1: For the examiner to establish the sequence of colors	4
Table 2: Data sheet	5
Ordering Information	6

DESCRIPTION:

For optometric, clinical and laboratory experimentation work. The Perimeter is used to map the areas of the retina sensitive to stimuli of different colors. Included with the apparatus are stimulus patches of red, yellow, green, blue, and white, 1/8" square, on a 5/8" aluminum disk.

The Perimeter is constructed in such a manner that measurements can be made of the threshold between regions where the stimulus is and is not discriminated. This is accomplished by having Subject (S) fixate on central spot (small mirror) while Experimenter (E) introduces the stimulus patch by moving it along the arm of the Perimeter. The distance of the stimulus from the eye is held constant by moving the stimulus target along the slit opening of the curved metal arm of the Perimeter, the eye being located at the center of the curvature of the arm. By moving the stimulus along the arm, E establishes the boundary between the regions of successful and unsuccessful color discrimination. The spatial position of the threshold is expressed by the angle between the line of fixation and the line from the eye to the threshold. Furthermore, by rotating the Perimeter arm about its pivot, threshold measurements along any line radiating out from the center of the visual field can be made; and by moving the stimulus patch inward and outward, threshold measurements for direction of stimulus movement can be made.

PROCEDURE:

Give S a comfortable seat and place cheekbone against back plastic support so that his left eye is directly centered in the fixation mirror. Have S keep his right eye covered during the experiment.

The nine rows of Table 2 show the angular positions of the Perimeter arm at which measurements are to be made. The columns give the hues of the stimulus patches and specify the direction in which E must move the stimulus.

In order that S be as naïve as possible with regard to the color being used, employ the following table or randomized stimulus sequences.

Table 1: A Table for the Experimenter to Establish the Sequence of Colors

Row	Sequence	Row	Sequence
1	G Y B W R	5	G Y R B W
2	B R W G Y	6	R B G W Y
3	Y B R W G	7	B W R Y G
4	Y W G B R	8	R G Y W B

Choose a row from Table 1—say Row 3—with the arm set at 0 degrees, present the stimuli in the order shown i.e., yellow, blue, red, white, and green. Present each color twice—first moving inward, then outward, beginning with the stimulus patch at the periphery and *slowly* moving the patch into the field of vision. The s must have no foreknowledge of the color for the *in* procedure. Tell S to report aloud and continuously what he sees—for example, “nothing, nothing, gray, blue, blue, yellow.” When the color is identified, stop and then read and record the angular distance from the center of fixation under *in* and *yellow* in the first row of Table 2. Next, administer the *out* trial, starting well within S’s field of color vision, and slowly moving the yellow patch toward the periphery. Again instruct S to make a continuous oral report. Do not stop movement until S reports a change in the color.

After completing the above procedure for the five colors, rotate the metal arm 45 degrees to obtain the data required for the second row of Table 2. The sequence of colors may again be determined by a random selection from Table 1 or you may simply proceed to the next available row returning to Row 1 following the completion of Row 8.

Table 2
Data Sheet

Angle of Rotation	Red		Yellow		Green		Blue		White	
	In	Out	In	Out	In	Out	In	Out	In	Out
0										
45										
90										
135										
180										
225										
270										
315										
360										

Lafayette Instrument Stand Perimeter Schweigger Type

Model 14014 User's Manual

Ordering Information:

All phone orders must be accompanied by a hard copy of your order. All must include the following information:

- 1) Complete billing and shipping addresses
- 2) Name and department of end user
- 3) Model number and description of desired item(s)
- 4) Quantity of each item desired
- 5) Purchase order number or method of payment
- 6) Telephone number

DOMESTIC TERMS

There is a \$50 minimum order. Open accounts can be extended to most recognized educational institutions, hospitals and government agencies. Net amount due 30 days from the date of shipment. Enclose payment with the order; charge with VISA, MasterCard; or pay COD. We must have a hard copy of your order by mail or fax. Students, individuals and private companies may call for a credit application.

INTERNATIONAL PAYMENT INFORMATION

There is a \$50 minimum order. Payment must be made in advance by: draft drawn on a major US bank; wire transfer to our account; charge with VISA, MasterCard; or confirmed irrevocable letter of credit. Proforma invoices will be provided upon request.

RETURNS

Equipment may not be returned without first receiving a Return Goods Authorization Number (RGA).

When returning equipment for service, please call Lafayette Instrument to receive a RGA number. Your RGA number will be good for 30 days. Address the shipment to: Lafayette Instrument Company, 3700 Sagamore Parkway North, Lafayette, IN 47904, U.S.A. Shipments cannot be received at the PO Box. The items should be packed well, insured for full

value, and returned along with a cover letter explaining the malfunction. Please also state the name of the Lafayette Instrument representative authorizing the return. An estimate of repair will be given prior to completion ONLY if requested in your enclosed cover letter. We must have a hard copy of your purchase order by mail or fax, or repair work cannot commence.

WARRANTY

Lafayette Instrument guarantees its equipment against all defects in materials and workmanship to the ORIGINAL PURCHASER for a period of one (1) year from the date of shipment, unless otherwise stated. During this period, Lafayette Instrument will repair or replace, at its option, any equipment found to be defective in materials or workmanship. If a problem arises, please contact our office for prior authorization before returning the item. This warranty does not extend to damaged equipment resulting from alteration, misuse, negligence or abuse, normal wear or accident. In no event shall Lafayette Instrument be liable for incidental or consequential damages. There are no implied warranties or merchantability of fitness for a particular use, or of any other nature. Warranty period for repairs or used equipment purchased from Lafayette Instrument is 90 days.

DAMAGED GOODS

Damaged equipment should not be returned to Lafayette Instrument prior to thorough inspection.

When a shipment arrives damaged, note damage on delivery bill and have the driver sign it to acknowledge the damage. Contact the delivery service, and they will file an insurance claim. When damage is not detected at the time of delivery, contact the carrier and request an inspection within 10 days of the original delivery. Please call the Lafayette Instrument Customer Service Department for a return authorization for repair or replacement of the damaged merchandise.



3700 Sagamore Parkway North
P.O. Box 5729 • Lafayette, IN 47903 USA
Tel: 765.423.1505 • 800.428.7545
Fax: 765.423.4111
E-mail: lic@lafayetteinstrument.com
www.lafayetteinstrument.com

Lafayette Instrument Co. Europe

4 Park Road, Sileby,
Loughborough, Leics., LE12 7TJ. UK.
Tel: +44 (0)1509 817700
Fax: +44 (0)1509 817701
E-mail: EUsales@lafayetteinstrument.com