

## Stand Perimeter, Schweiger Type Model 14014



### Product Description

Subject fixes eyes at midpoint in the perimeter's semicircular while experimenter moves a pointer to identified by subject at his / her visual periphery. This device is used to map areas of the retina sensitive to different colored stimuli. Stand perimeter comes with semicircular perimeter device with fixation point, five targets (blue, green, yellow, red, and white), a wand for presenting targets, five right eye charts, five left eye charts, and stand.

### Optional Products for Usage

#### Set of 5 Extra Targets

Model 14100

Replacement Set of Targets for the Stand Perimeter, Model 14014.

#### Replacement Wand

Model 14101

Replacement Wand for the Stand Perimeter, Model 14014.

#### Right Eye Charts

Model 14102

Replacement Right Eye Charts for the Stand Perimeter, Model 14014.

#### Left Eye Charts

Model 14103

Replacement Left Eye Charts for the Stand Perimeter, Model 14014.

---

### Contact Information

S.I. Instruments  
256 South Rd. Hilton  
South Australia 5033  
Ph (08) 8352 5511

[info@si-instruments.com](mailto:info@si-instruments.com)  
[www.si-instruments.com](http://www.si-instruments.com)