

S.I. Instruments
256 South Rd. Hilton
South Australia 5033
Ph (08) 8352 5511

info@si-instruments.com.au
si-instruments.com.au

Manual Dexterity

Rehabilitation and Return to Work

Skills Assessment

Perception

Motor Coordination

Personnel Selection and Development



Purdue Pegboard Test

Model 32020A

The Purdue Pegboard Test was first developed by Joseph Tiffin, Ph.D., an Industrial Psychologist at Purdue University in 1948. Since that time, this device has been used extensively to aid in the selection of employees for jobs that require fine and gross motor dexterity and coordination. It measures gross movements of hands, fingers, arms, and fingertip dexterity as necessary in assembly tasks.

Intended for industrial use and assembly work in a factory setting, the test is now being studied for use in other special areas such as patients with Parkinson's, Multiple Sclerosis, Stroke sufferers, and similar illnesses. Physical and Occupational Therapists also use the Purdue Pegboard for injury rehabilitation, using the test as a tool to obtain baseline data on a patient documenting the progress and/or degree of disability.


Complete with the following tools

The Pegboard is equipped with pins, collars, and washers which are located in four cups at the top of the board. An examiner's manual with norms is also included for administering the test.

Replacement Models

32103 Replacement set of 55 pegs, 45 washers, and 25 collars

32107 Record blanks, 25 per package



Scan for more
information on the
Purdue Pegboard

Purdue Pegboard Scoring App

Lafayette Instrument is proud to unveil the Purdue Pegboard Scoring Application for Android™ and iOS. This tablet-based app is intended for use with the Lafayette Instrument Company Purdue Pegboard Test (Models 32020 and 32020A), and it assists administrators in all areas of the testing process by standardizing administration through easily set up test batteries with optionally read instructions, creating organizational norms, and keeping track of individualized data.



**Visit Google Play or the iOS App Store
today to download the free 12 test trial.***

* In app purchases required for extended usage.

Grooved Pegboard Test

Model 32025


This manipulative dexterity test contains twenty-five holes with randomly positioned slots and pegs which have a key along one side. Pegs must be rotated to match the hole before they can be inserted.

This procedure measures performance speed in a fine motor task by examining both sides of the body, inferences may be drawn regarding possible lateral brain damage.

The test requires more complete visual motor coordination than most of our pegboards and has been used in several neuropsychological test batteries, in student labs, and as a screening technique in industrial environment.

The Grooved Pegboard is equipped with pegs, and examiner's manual with norms.

Replacement Model: 32104 Replacement Pegs, 32 per package



Hand Tool Dexterity Test

Model 32521

This test measures proficiency in using ordinary mechanical tools. The test consists of tools and two uprights with bolts, washers, and nuts. The object is to disassemble all the bolts from one upright and reassemble them on corresponding rows of the other upright with the heads of the bolts inside.

This type of skill is important to many industrial jobs and apprentice training. Results of the test have been used to determine vocational interest and as an indicator of success where job/tasks require the use of these or similar tools.

The Hand Tool Dexterity Test is complete with examiner's manual, norms, 15/16 inch Open End Wrench, 1/2 inch Open End Wrench, 10 inch Crescent Wrench, Screw Driver, and all required Bolts, Nuts, and Washers.



Replacement Models

- 32521A 15/16 inch Open End Wrench
- 32521B 1/2 inch Open End Wrench
- 32521C 10 inch Crescent Wrench
- 32521D Screwdriver
- 32521E Bolts, Nuts, and Washers

Minnesota Manual Dexterity Tests

Minnesota Manual Dexterity Test Model 32023

The Minnesota Manual Dexterity testing kit includes sixty plastic cylinders, 1 board with sixty round holes, carrying case, manual with norms, and record blanks.

The Minnesota Manual Dexterity Test consists of the following two tests:

- Placing
- Turning

Complete Minnesota Manual Dexterity Test Model 32023A

The Complete Minnesota Manual Dexterity Test kit includes 2 boards, 60 blocks, carrying case, and instruction manual with norms.

The Complete Minnesota Manual Dexterity Test consists of a battery of five tests:

- Placing
- Turning
- Displacing
- One-Hand Turning and Placing
- Two-Hand Turning and Placing

This widely used test measures capacity for simple but rapid hand-eye coordination. This is particularly applicable in shop and office occupations requiring quick movement in handling simple tools and production materials without differentiating size and shape.

Replacement Models

1-04811 Replacement Board

32031SET Replacement cylinders (complete set)

32031 Replacement cylinders (sold individually)

32032 Record blanks for 32023, 50 per package

32033 Record blanks for 32023A, 50 per package



O'Connor Finger Dexterity Test

Model 32021

The O'Connor Finger Dexterity Test requires hand placement of 3 pins per hole. Consisting of 100 3/16" diameter holes with holes arranged in ten rows and spaced ½ inch apart. Primarily used as a predictive tool wherever rapid manipulation of objects, especially the picking up and placing of small parts, is important.

The O'Connor Finger Dexterity Test is complete with board, pins, examiner's manual, and norms.

Replacement Model: 32106 Set of Replacement Pins, 105 per package



O'Connor Tweezer Dexterity Test

Model 32022

This test requires the use of tweezers to place a single pin in each 1/16" diameter hole. A high score is indicative of manual aptitude for work involving precision and steadiness of small hand tools and requires a high degree of hand-eye coordination.

The O'Connor Tweezer Dexterity Test is complete with board, pins, tweezers, examiner's manual, and norms.

Replacement Models

32106 Set of Replacement Pins, 105 per package

32109 Replacement Tweezers



Two Arm Coordination Test

Model 32532

This is a test of the coordination of both arms working together in order to move a stylus around a six-point star pattern. This unit must be connected to an impulse counter to record the number of errors and/or a stop clock to record the amount of time outside the path. Norms not included.

Available Counters

54060A Clock/Counter is a timing device that is capable of accurately recording times down to a millisecond.



Roeder Manipulative Aptitude Test

Model 32026

This test measures hand, arm, finger dexterity, and speed. It is designed to test individuals for employment and to test elementary through college students when dexterity is a primary requirement.


The board has four receptacles for holding washers, rods, caps, and nuts. The performance board also is comprised of a horizontal T-bar and 40 inserts arranged in a predetermined pattern.

The Roeder Manipulative Aptitude Test is complete with examiner's manual, norms, board, as well as all washers, rods, caps, and nuts.

Replacement Models

32026RB Set of Record Blanks, 50 per package

32026P Replacement Washers, Rods, Caps, and Nuts



Auto Scoring Mirror Tracer

Model 58024E

This tracing device involves reversal ability, hand-eye coordination, and learning. The subject is required to trace the star pattern while watching only its mirror image. Errors are automatically tallied and provision has also been made to track the total amount of time that the stylus is outside the star pattern. The mirror and shield are completely collapsible for easy storage and portability. The system has been designed to be used with either left or right handed subjects without a need to change the mounting brackets, support arms, etc.

This unit must be connected to a timing device to also record the amount of time outside of the path. Norms are not included.

Available Timers


54060A Clock/Counter is a timing device that is capable of accurately recording times down to a millisecond.

Replacement Parts

32532A Replacement Star

32533B Replacement Stylus

58024M Replacement Mirror and Base




Occupational Skills Assessment Test Battery

Model 32604

This test battery is used to measure the progress in rehabilitation and return-to-work capability of individuals performing jobs/tasks that require manual dexterity, hand-eye coordination, steadiness, and perceptual motor skills. Placement personnel and human resource departments can pretest job applicant suitability for assembly and other jobs where frequent manipulation of objects in confined spaces is required.

Complete with the following tools (photos not to scale)


- | | |
|---------------------------------------|---|
| 32011 Steadiness Tester, Hole Type | 32026 Roeder Manipulative Aptitude Test |
| 32020A Purdue Pegboard Test | 32521 Hand Tool Dexterity Test |
| 32022 O'Connor Tweezer Test | 32532 Two-Arm Coordination Test |
| 32023 Minnesota Manual Dexterity Test | 58024E Auto-Scoring Mirror Tracer |



32020A



32011



32026



32521




32532



32022



32023



58024E



3700 Sagamore Pkwy N
Lafayette, IN 47904

Phone: (765) 423-1505
Fax: (765) 423-4111

sales@lafayetteinstrument.com
www.lafayetteevaluation.com



S.I. Instruments
256 South Rd. Hilton
South Australia 5033
Ph (08) 8352 5511

info@si-instruments.com.au
si-instruments.com.au